




**BIOCREDIT**

One Step

**COVID-19 IgG+IgM Duo**



ANTIBODY DETECTION TEST  
KIT

# BIOCREDIT COVID-19 IgM / IgG DUO

명칭(제품명, 품목명, 모델명): BIOCREDIT COVID-19 IgG+IgM Duo, 고위험서감염  
포장단위: 25 tests/kit  
이허용, 사용 금지

- Storage: Store at 1~40°C (33.8~104°F)
- Package: 25 tests/kit
- Do not reuse

LOT Lot No: HD 7  
Exp. 20  
REF G66 A25






# VIRUS AND VIRUS DIAGNOSTICS

According to the statistics of **Rospotrebnadzor** almost **30%** of those infected with coronavirus in Russia **are asymptomatic**.

**Mass testing** of the population for COVID-19 which is recommended by the World Health Organization is a proven and effective infection-control tool.






# MEDICAL DEVICE PURPOSE

Designed for qualitative detection of IgM and IgG antibodies to COVID-19 coronavirus in the whole blood, serum or plasma samples.


When using the Kit for biological samples testing the results indicating the sample reactivity should be confirmed by alternative test methods.

The kit is suitable for in vitro use only.



# ANTIBODIES TO COVID-19

Graph of the antibodies detection to SARS-CoV-2 with relation to the onset of the disease.



## IgM ANTIBODIES

indicate the recent infection and can be used as an auxiliary diagnosis of the early infection.

## IgG ANTIBODIES

are produced later and for a longer period and indicate of a previous or a secondary infection.

Source: [gemohelp.ru](http://gemohelp.ru)

# INDICATIONS

For the qualitative determination of the presence of IgM or IgG antibodies to coronavirus (SARS-CoV-2) in the serum, plasma and whole blood samples in patients with clinical symptoms of respiratory disease with suspected COVID-19 infection, as well as in patients without the signs of colds and not in contact with the COVID-19 patients.

# CONTRAINDICATIONS

- Expired test.
- Damaged container with the device.
- Inadequate product storage and transportation.
- There are no other contraindications unless the specimen cannot be taken for medical reasons.

# TEST KIT FEATURES

- Run time - 5-10 minutes;
- Auxiliary tool for diagnosis;
- Easy-to-use, test strips (cassettes).
- Identification of those who were ill and acquired immunity.





# KIT CONTENTS

BIOCREDIT COVID-19 IgG/IgM DUO kit contains:



- Test cassette (in individual foil package with desiccant material) – 25 pieces (IgG) + 25 pieces (IgM).
- A vial with the solution for biological sample diluting — 1 pc.
- Disposable capillary pipettes – 50 pc.
- Prescribing information - 1 pc.

Materials not included in the scope of supply: protective gloves, protective glasses, micropipette, alcohol soaked wipe, scarifier(lancet), blood sampling tubes.

# CLINICAL AND LABORATORY TESTING IGm

Marketing Authorization issued by Roszdravnadzor  
Ref. No.2020/10777 dated June 10, 2020.

93.4%

TEST SENSITIVITY

100%

DIAGNOSTIC  
SPECIFICITY






EDIT COVID-19 IgG+IgM Duo

Marketing Authorization issued by Roszdravnadzor  
Ref. No.2020/10777 dated June 10, 2020.

# CLINICAL AND LABORATORY TESTING IgG




TEST SENSITIVITY



DIAGNOSTIC  
SPECIFICITY





## TESTING PROCEDURE SCHEME

1. Add 10  $\mu\text{l}$ (serum or plasma) or 20  $\mu\text{l}$  (whole blood) of the sample to sample well (S), using the capillary pipette or micropipette included in the set.
2. Add 3 drops of the diluent for analyses to buffer well (B) of the device.
3. Read the result during 5~10 minutes.




# INTERPRETATION OF RESULTS

If only the control line (C) is colored red in the result window: the test is negative for the presence of the coronavirus antigen.


If the control line (C) is colored red and the test line (T) turns black in the result window: the test is positive for the presence of the coronavirus antigen.

If there is no red-colored control line (C) in the result window, the test result is considered invalid. The reason for such a result could be a violation of the biological sample collection and/or analysis procedures or the failure of the test cassette(test system) used.


Note: strip width or color intensity are irrelevant.



Positive



Negative



Invalid





# OUR CLIENTS



**ПРОДО**

