REAGENT KIT FOR
IMMUNOCHROMATOGRAPHIC TEST FOR
SARS-COV-2 ANTIGEN DETECTION IN
NASOPHARYNGEAL SWABS
(ND COVID-19 AG TEST).

MEDICAL DEVICE PURPOSE

Designed for the qualitative detection of SARS-CoV-2 antigen in human biological material (nasopharyngeal swab).

The kit is an auxiliary tool for early SARS CoV-2 diagnosis in patients with clinical symptoms of SARS-CoV-2 infection. When biological samples are tested using the Kit, the results indicating that the sample is reactive should be confirmed by alternative methods of analysis.

The kit is intended for laboratory in vitro diagnostics only.

INDICATIONS

Designed for the qualitative detection of the SARS-CoV-2 antigen in nasopharyngeal swab in patients with respiratory disease symptoms with a suspected COVID-19 infection.

CONTRAINDICATIONS

- Expired test
- Damaged package
- Inadequate product storage and transportation.
- There are no other contraindications unless the specimen cannot be taken for medical reasons.

PRINCIPLES OF THE TEST

The reagent kit action is based on the principle of using a two-color system. Medical device structure includes a substrate wherein colloidal gold particles are included to form conjugated complexes, and a membrane strip of reagents where antibodies specific to the SARS-CoV2 coronavirus antigen are preliminarily applied.

The test strip (T) is coated with mouse monoclonal antibodies to the SARS-CoV-2 coronavirus antigen.


The control line (C) is coated with monoclonal rabbit antibodies to chicken yolk immunoglobulin IgY. The SARS-CoV-2 antigen in a biological sample during the test reacts with the mouse anti-SARS-CoV-2 monoclonal antibody conjugated with dye particles to form an antigen-antibody colored particle complex. The complex reaches the test strip (T), where it is bound with mouse anti-SARS-CoV-2 monoclonal antibodies.

If the biological sample contains SARS-CoV-2 antigens, the color test line becomes visible with a naked eye in the result aperture.

TEST KIT FEATURES

- The total testing time does not exceed 20 minutes.
- Easy-to-use test strips (cassettes)

KIT CONTENTS

The device is presented by the set of reagents, which includes:

1. Test cassette (together with a desiccant in individual foil package) - 25 pcs.;
2. Test tube with biological sample extraction buffer - 25 pcs.;
3. Cardboard stand - 1 pc.;
4. Prescribing information - 1 pc.

The following materials are not provided

- Protective gloves;
- Timer or stopwatch;
- Micropipettes;
- Sterile swab for specimen collection

CLINICAL AND LABORATORY TESTING




96%

TEST SENSITIVITY

98%

DIAGNOSTIC SPECIFICITY

INTERPRETATION OF RESULTS

Note: strip width and color intensity are irrelevant.

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