

**REAGENT KIT FOR
IMMUNOCHROMATOGRAPHIC TEST FOR
SARS-COV-2 ANTIGEN DETECTION**

BIOCREDIT COVID-19 Ag

MEDICAL DEVICE PURPOSE

Designed for qualitative detection of SARS-CoV-2 antigen in a nasopharyngeal swab.

The kit is suitable for in vitro use only.



INDICATIONS

For the qualitative detection of the SARS-CoV-2 antigen in nasopharyngeal swab specimens from individuals with respiratory disease signs and symptoms who are suspected of COVID-19.

CONTRAINDICATIONS

- Expired test.
- Damaged container with the device.
- Inadequate product storage and transportation.
- There are no other contraindications unless the specimen cannot be taken for medical reasons.

REAGENT KIT FEATURES

BIOCREDIT COVID-19 Ag Reagent Kit is designed for the qualitative and rapid detection of the COVID-19 antigen in nasopharyngeal swab.

Easy-to-use test strips.

Testing time 5 minutes.



CLINICAL AND LABORATORY TESTING

96%

TEST SENSITIVITY

100%

DIAGNOSTIC SPECIFICITY



PRINCIPLES OF THE TEST

The reagents kit for SARS-CoV2 antigen detection (BIOCREDIT COVID-19 Ag) is designed for **rapid and quality** immunochromatographic assay.

T LINE

Antibodies specific to the SARS-CoV-2 coronavirus antigen were preliminarily applied to the membrane strip of reagent **T Line**. If the test biological sample contains the SARS-CoV-2 coronavirus antigen, a visible black strip appears along the line at the result of testing.

C LINE

C Control line is provided in the test system to provide integration of the correct test execution procedure. If the test procedure is performed correctly, control line C should be colored under all circumstances, regardless of the test results.

KIT CONTENTS




The product is a reagent kit containing:

- Test cassette (together with a desiccant in individual foil package) - 20 pcs.
- A test tube with a buffer for biological sample diluting - 20 pcs.
- Dropper cap - 20 pcs.
- Prescribing information - 1 pc.

The following materials are not provided: protective gloves, safety goggles.

SAMPLE COLLECTION



Biological samples should be taken by trained healthcare professionals. Biological samples should be handled with caution due to the risk of infection.

Handling of biological samples requires special care; any violation of prescribed procedures may adversely affect the specimen and lead to inaccurate test results.


In order to obtain more accurate test results, it is advisable to take a biological sample from several areas of the nasopharyngeal.

IMPORTANT

Test should be performed immediately as soon as a sample is taken. If this is not possible, a swab with the sampled material may be stored for a maximum 12 hours in a refrigerator, at a temperature within 2 - 8°C, or for a maximum of 24 hours in a freezer at -20°C.



NASOPHARYNGEAL SWAB COLLECTION



To collect a biological sample from the nasopharyngeal area, carefully insert a swab on a stick into the nasal cavity.

Keep gently inserting the tampon until resistance is encountered in the nasal turbinate.

Gently rotate the swab around its axis several times with your fingers and remove it from the nasal cavity. Visually verify that the tip of the swab is wet.

TEST PROCEDURE

PREPARATORY STEP

Ensure that the biological sample and all Kit components are at room temperature before performing the test.

Do not open an individual test cassette packaging until all preliminary manipulations are completed and test procedure is ready for immediate performance.

TEST STAGES

Peel off aluminum foil seal from the top of the extraction vial containing the extraction buffer.

Place a swab with a smear from the nasopharyngeal area into the extraction vial with a buffer and press the swab tip against inner wall or bottom of the vial to release biological sample and transfer it to the buffer medium.

While continuing to press the swab tip against the wall or bottom, make 5 – 10 circular movements.

Press the tip of the swab against the inner wall to release the liquid, then remove the swab from the vial and dispose it in accordance with the hazardous biological waste disposal requirements.

Securely close the vial using the dropper cap from the kit.

Remove the test cassette from the individual package. Place the cassette on a dry, flat surface.



TEST STAGES

Turn down the vial with the buffer and biological sample and press lightly on the sides to dispense 3 - 4 drops (90 - 150 μ l). Obtained volume should be placed in the center of the sample well on the surface of the test cassette, marked with the letter S. Ensure that an adequate volume of buffer containing biological sample is used for the test.

Insufficient or excessive sample volume, as well as insufficient or excessive buffer volume, can lead to inaccurate (invalid or distorted) test results.

Read the test result in 5 - 8 minutes.






TESTING PROCEDURE SCHEME


Read the test result in 5 - 8 minutes.

INTERPRETATION OF RESULTS


Note: strip width or color intensity are irrelevant.




If only the control line (C) is colored red in the result window: the test is negative for the presence of the coronavirus antigen.



If the control line (C) is colored red and the test line (T) turns black in the result window: the test is positive for the presence of the coronavirus antigen.



If there is no red-colored control line (C) in the result window, the test result is considered invalid. The reason for such a result could be a violation of the biological sample collection and/or analysis procedures or the failure of the test cassette(test system) used.



LIMITATIONS

The test result can be negative if the coronavirus content (titer) in the biological sample is below the detection limit.


A negative test result does not rule out the presence of infection in cases of recent infection.



CERTIFICATES

DECLARATION OF CONFORMITY

Manufacturer: RapiGEN Inc.
2F, 25, Heungan-daero, Gunpo-si, Gyeonggi-do 15809, Republic of Korea

European Representative: MT Promed Consulting GmbH
Altenhofstrasse, 80, 66586 St. Ingbert, Germany 

Product: BIOCREREDIT COVID-19 Ag
Catalog no.: G61RHA20

Classification: Neither listed in the annex II of the IVDD, nor self-testing device
EDMA code: 15.70.90.90.00; Other Other Virology Rapid Tests

Conformity Assessment Route: Self Declaration (according to annex III of IVDD)

We herewith declare that the above mentioned products meet the provisions of the council Directive 93/79EC for in vitro diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer.

Standards applied: ISO13485:2016, EN ISO14971:2012, EN13640:2002, EN13641:2002, EN13612:2002, EN ISO 15223-1:2016, ISO17511:2003, EN13975:2003, EN ISO 18113-1:2011, EN ISO 18113-2:2011, EN 62366-2008, KGMP

Place, Date of Issue: Gyeonggi-Do, Republic of Korea, 01st April, 2020.

Signature: 
Jae-Ku, Park
CEO/President
RapiGEN Inc.



Jae-Ku Park, President





PUBLICATIONS IN MEDIA

MEDRUSSIA.ORG portal- a rapid test for SARS-CoV-2 antigen detection has appeared in Russia: Biocredit from RapiGEN.

GOVORIT MOSKVA radio station- Moscow ambulances will be equipped with an innovative Korean test for coronavirus detection

ARGUMENTY I FAKTY- Yelena Malysheva: “Know your coronavirus test in 5 minutes”

VADEMECUM - THE MINISTRY OF HEALTH RECOMMENDED TO USE EXPRESS TESTS FOR COVID-19 "AT THE PATIENT'S SIDE" AND MORE

MEDVESTNIK.RU portal- according to Korean system

VADEMECUM - SOUTH KOREAN RAPIGEN EXPRESS TEST FOR COVID-19 OBTAINED MARKETING AUTHORIZATION IN RUSSIA